


Galectin 3 Antibody (R33294)

Catalog No.	Formulation	Size
R33294-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3958
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This Galectin 3 antibody is available for research use only.



Western blot testing of human [A] and mouse [B] colon lysates with Galectin 3 antibody at 0.3ug/ml. Predicted molecular weight ~26kDa.

Description

Additional name(s) for this target protein: Lectin, galactoside-binding, soluble, 3; LGALS3

Application Notes

Optimal dilution of the Galectin 3 antibody should be determined by the researcher.

Immunogen

Amino acids CNTKLDNNWGREERQ were used as the immunogen for this Galectin 3 antibody.

Storage

Aliquot and store the Galectin 3 antibody at -20oC.