


## Galectin-3 Antibody (R34764)

| Catalog No.  | Formulation                                                                   | Size   |
|--------------|-------------------------------------------------------------------------------|--------|
| R34764-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

|                    |                                                              |
|--------------------|--------------------------------------------------------------|
| Availability       | 1-3 business days                                            |
| Species Reactivity | Mouse                                                        |
| Format             | Antigen affinity purified                                    |
| Clonality          | Polyclonal (goat origin)                                     |
| Isotype            | Goat Ig                                                      |
| Purity             | Antigen affinity                                             |
| Gene ID            | 16854                                                        |
| Applications       | Western Blot : 0.5-2ug/ml<br>ELISA (peptide) LOD : 1:4000    |
| Limitations        | This Galectin-3 antibody is available for research use only. |



Western blot testing of mouse heart lysate with Galectin-3 antibody at 0.5ug/ml.  
Predicted molecular weight ~26kDa.

## Description

Additional name(s) for this target protein: Lectin, galactose binding, soluble 3; Lgals3

## Application Notes

Optimal dilution of the Galectin-3 antibody should be determined by the researcher.

## Immunogen

Amino acids NANRIVLDFRRGND were used as the immunogen for this Galectin-3 antibody.

## Storage

Aliquot and store the Galectin-3 antibody at -20oC.