


## GATA6 Antibody (F49005)

| Catalog No.   | Formulation                                | Size    |
|---------------|--------------------------------------------|---------|
| F49005-0.4ML  | In 1X PBS, pH 7.4, with 0.09% sodium azide | 0.4 ml  |
| F49005-0.08ML | In 1X PBS, pH 7.4, with 0.09% sodium azide | 0.08 ml |

[Bulk quote request](#)

|                           |                                                         |
|---------------------------|---------------------------------------------------------|
| <b>Availability</b>       | 1-3 business days                                       |
| <b>Species Reactivity</b> | Human                                                   |
| <b>Format</b>             | Antigen affinity purified                               |
| <b>Clonality</b>          | Polyclonal (rabbit origin)                              |
| <b>Isotype</b>            | Rabbit Ig                                               |
| <b>Purity</b>             | Antigen affinity                                        |
| <b>UniProt</b>            | Q92908                                                  |
| <b>Applications</b>       | Western Blot : 1:1000                                   |
| <b>Limitations</b>        | This GATA6 antibody is available for research use only. |



### Description

GATA6 is thought to be important for regulating terminal differentiation and/or proliferation.

### Application Notes

Titration of the GATA6 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

### Immunogen

A portion of amino acids 254-279 from the human protein was used as the immunogen for this GATA6 antibody.

### Storage

Aliquot the GATA6 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.