


GBA Antibody / Glucosylceramidase beta / GC [clone AACCA-7] (RQ5035)

Catalog No.	Formulation	Size
RQ5035	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul

[Bulk quote request](#)

Availability	1-2 weeks
Species Reactivity	Human
Format	Purified
Clonality	Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	AACA-7
Purity	Affinity purified
UniProt	P04062
Applications	Western Blot : 1:500-1:2000
Limitations	This GBA antibody is available for research use only.



Western blot testing of human U-87 MG cell lysate with GBA antibody. Predicted molecular weight ~60 kDa.

Description

This gene encodes a lysosomal membrane protein that cleaves the beta-glucosidic linkage of glycosylceramide, an intermediate in glycolipid metabolism. [RefSeq]

Application Notes

Optimal dilution of the GBA antibody should be determined by the researcher.

Immunogen

A synthetic peptide specific to human Glucosylceramidase / GBA was used as the immunogen for the GBA antibody.

Storage

Store the GBA antibody at -20oC.