


GCNT3 Antibody (R34837)


Catalog No.	Formulation	Size
R34837-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	9245
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This GCNT3 antibody is available for research use only.



Western blot testing of human duodenum lysate with GCNT3 antibody at 0.1ug/ml. Predicted molecular weight ~51kDa.



Western blot of HEK293 lysate overexpressing GCNT3 probed with GCNT3 antibody (mock transfection in lane 1). Predicted molecular weight ~51kDa.

Description

Additional name(s) for this target protein: Glucosaminyl (N-acetyl) transferase 3, C2GnT-M

Application Notes

Optimal dilution of the GCNT3 antibody should be determined by the researcher.

Immunogen

Amino acids NKFDPKVDDNALQ were used as the immunogen for this GCNT3 antibody.

Storage

Aliquot and store the GCNT3 antibody at -20oC.