

Gelsolin Antibody (R35115)

Catalog No.	Formulation	Size
R35115-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2934
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This Gelsolin antibody is available for research use only.



Western blot testing of human duodenum lysate with Gelsolin antibody at 1ug/ml.
Predicted molecular weight ~86kDa, observed here at ~100kDa.

Description

Additional name(s) for this target protein: GSN, ADF, AGEL

Application Notes

Optimal dilution of the Gelsolin antibody should be determined by the researcher.

1. This Gelsolin antibody is specific for isoform a.

Immunogen

Amino acids PQGRVPEARPNS were used as the immunogen for this Gelsolin antibody.

Storage

Aliquot and store the Gelsolin antibody at -20oC.