


GGT3P Antibody (Center Region) (F54064)

Catalog No.	Formulation	Size
F54064-0.2ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.2 ml
F54064-0.05ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.05 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	A6NGU5
Applications	Western Blot : 1:2000
Limitations	This GGT3P antibody is available for research use only.



Description

Initiates extracellular glutathione (GSH) breakdown; catalyzes the transfer of the glutamyl moiety of glutathione to amino acids and dipeptide acceptors. [UniProt]

Application Notes

The stated application concentrations are suggested starting points. Titration of the GGT3P antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 176-207 from human GGT3P was used as the immunogen for the GGT3P antibody.

Storage

Aliquot the GGT3P antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.