


GIRK2 Antibody / KCNJ6 (R36067)

Catalog No.	Formulation	Size
R36067-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3763
Applications	Western Blot : 0.5-2ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This GIRK2 antibody is available for research use only.



Western blot testing of human brain (substantia nigra) lysate with GIRK2 antibody at 0.5ug/ml. Predicted molecular weight ~48 kDa.

Description

Additional name(s) for this target protein: G protein-activated inward rectifier potassium channel 2; KCNJ6; potassium inwardly-rectifying channel, subfamily J

Application Notes

Optimal dilution of the GIRK2 antibody should be determined by the researcher.

Immunogen

Amino acids SSKLNQHAELET were used as the immunogen for this GIRK2 antibody.

Storage

Aliquot and store the GIRK2 antibody at -20oC.