


Glypican-1 Antibody / GPC1 (R36436)

Catalog No.	Formulation	Size
R36436-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2817
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This Glypican-1 antibody is available for research use only.



Western blot testing of human 1) HeLa and 2) K562 cell lysate with Glypican-1 antibody at 2ug/ml. Predicted molecular weight: ~62 kDa.

Description

Additional name(s) for this target protein: Glypican proteoglycan 1

Application Notes

Optimal dilution of the Glypican-1 antibody should be determined by the researcher.

Immunogen

Amino acids TDKFWGTSGVES were used as the immunogen for this Glypican-1 antibody.

Storage

Aliquot and store the Glypican-1 antibody at -20oC.