


GNAT3 Antibody (R34552)

Catalog No.	Formulation	Size
R34552-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	346562
Applications	Immunohistochemistry (FFPE) : 3-6ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This GNAT3 antibody is available for research use only.



IHC testing of FFPE human small intestine with GNAT3 antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining. Staining of endocrine cells amongst the enterocytes is seen.

Description

Additional name(s) for this target protein: Guanine nucleotide binding protein, alpha transducing 3; GNAT3

Application Notes

Optimal dilution of the GNAT3 antibody should be determined by the researcher.

Immunogen

Amino acids KNQFLDLNLKKEDKE were used as the immunogen for this GNAT3 antibody.

Storage

Aliquot and store the GNAT3 antibody at -20oC.