


GNRHR Antibody (R35341)

Catalog No.	Formulation	Size
R35341-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2798
Applications	Western Blot : 0.5-2ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This GNRHR antibody is available for research use only.



Western blot testing of human peripheral blood mononucleocyte lysate with GNRHR antibody at 0.5ug/ml. Expected molecular weight: 38-60 kDa but may be observed at higher molecular weights due to glycosylation.

Description

Additional name(s) for this target protein: Gonadotropin-releasing hormone receptor, GRHR, LRHR

Application Notes

Optimal dilution of the GNRHR antibody should be determined by the researcher.

Immunogen

Amino acids QKWTQKKEKGK were used as the immunogen for this GNRHR antibody.

Storage

Aliquot and store the GNRHR antibody at -20oC.