


## GPBAR1 Antibody / TGR5 (RQ6095)

Catalog No.	Formulation	Size
RQ6095	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

**Bulk quote request**

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Mouse, Rat
<b>Format</b>	Antigen affinity purified
<b>Clonality</b>	Polyclonal (rabbit origin)
<b>Isotype</b>	Rabbit IgG
<b>Purity</b>	Affinity purified
<b>Buffer</b>	Lyophilized from 1X PBS with 2% Trehalose and 0.025% sodium azide
<b>UniProt</b>	Q8TDU6
<b>Applications</b>	Western Blot : 1-2ug/ml
<b>Limitations</b>	This GPBAR1 antibody is available for research use only.



Western blot testing of 1) rat liver, 2) mouse liver and 3) mouse kidney lysate with GPBAR1 antibody. Predicted molecular weight: 35-50 kDa depending on glycosylation level.

## Description

The G protein-coupled bile acid receptor 1 (GPBAR1) also known G-protein coupled receptor 19 (GPCR19), membrane-type receptor for bile acids (M-BAR) or TGR5 as is a protein that in humans is encoded by the GPBAR1 gene. This gene encodes a member of the G protein-coupled receptor (GPCR) superfamily. This enzyme functions as a cell surface receptor for bile acids. Treatment of cells expressing this GPCR with bile acids induces the production of intracellular cAMP, activation of a MAP kinase signaling pathway, and internalization of the receptor. The receptor is implicated in the suppression of macrophage functions and regulation of energy homeostasis by bile acids. Alternative splicing results in multiple transcript variants encoding the same protein.

## Application Notes

Optimal dilution of the GPBAR1 antibody should be determined by the researcher.

## Immunogen

Amino acids DQRYTAPWRAAAQR from the human protein were used as the immunogen for the GPBAR1 antibody.

## Storage

After reconstitution, the GPBAR1 antibody can be stored for up to one month at 4oC. For long-term, aliquot and store at -20oC. Avoid repeated freezing and thawing.