


GPR49 Antibody [clone 349CT16.1.3.1] (F40329)

Catalog No.	Formulation	Size
F40329-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F40329-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Clonality	Monoclonal (mouse origin)
Isotype	Mouse IgG1
Clone Name	349CT16.1.3.1
Purity	Purified
UniProt	O75473
Applications	Western Blot : 1?2000
Limitations	This GPR49 antibody is available for research use only.



Western blot analysis of K562 lysate and GPR49 antibody at 1:2000.

Description

Orphan receptor. Stem cell marker of the intestinal epithelium and the hair follicle. Target gene of Wnt signaling.

Application Notes

Titration of the GPR49 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 689-719 from the human protein was used as the immunogen for this GPR49 antibody.

Storage

Aliquot the GPR49 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.