


GPR83 Antibody (R35756)

Catalog No.	Formulation	Size
R35756-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Cow, Dog, Pig
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10888
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This GPR83 antibody is available for research use only.



Western blot testing of rat brain lysate with GPR83 antibody at 1ug/ml. Predicted molecular weight: ~48kDa.

Description

Additional name(s) for this target protein: G protein-coupled receptor 83

Application Notes

Optimal dilution of the GPR83 antibody should be determined by the researcher.

Immunogen

Amino acids RRYGAESQNPTVK were used as the immunogen for this GPR83 antibody. The amino acid sequence used

as immunogen is from the extracellular portion of the protein.

Storage

Aliquot and store the GPR83 antibody at -20oC.