


Granzyme B Antibody [clone RM441] (R20456)

Catalog No.	Formulation	Size
R20456-0.1ML	Antibody in PBS with 50% glycerol, 1% BSA and 0.09% sodium azide	100 ul


Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	RM441
Purity	Protein A purified from animal origin-free supernatant
UniProt	P10144
Localization	Cytoplasmic
Applications	Immunohistochemistry (FFPE) : 1:100-1:200
Limitations	This recombinant Granzyme B antibody is available for research use only.



IHC staining of FFPE human lung squamous carcinoma tissue with recombinant Granzyme B antibody at 1:200.



IHC staining of FFPE human tonsil tissue with recombinant Granzyme B antibody at 1:200.

Description

This antibody reacts to Granzyme B.

Application Notes

The stated application concentrations are suggested starting points. Titration of the recombinant Granzyme B antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A peptide corresponding to the C-terminus of Granzyme B was used as the immunogen for the recombinant Granzyme B antibody.

Storage

Store the recombinant Granzyme B antibody at -20oC.