


## GRIN3B Antibody (R33711)

| Catalog No.  | Formulation                                                                   | Size   |
|--------------|-------------------------------------------------------------------------------|--------|
| R33711-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

|                    |                                                             |
|--------------------|-------------------------------------------------------------|
| Availability       | 1-3 business days                                           |
| Species Reactivity | Human                                                       |
| Format             | Antigen affinity purified                                   |
| Clonality          | Polyclonal (goat origin)                                    |
| Isotype            | Goat Ig                                                     |
| Purity             | Antigen affinity                                            |
| Gene ID            | 116444                                                      |
| Applications       | Western Blot : 0.5-1.5ug/ml<br>ELISA (peptide) LOD : 1:8000 |
| Limitations        | This GRIN3B antibody is available for research use only.    |



Western blot testing of human cerebellum lysate with GRIN3B antibody at 0.5ug/ml.  
Predicted molecular weight: ~113kDa.

## Description

Additional name(s) for this target protein: Glutamate receptor, ionotropic, N-methyl-D-aspartate 3B

## Application Notes

Optimal dilution of the GRIN3B antibody should be determined by the researcher.

## Immunogen

Amino acids EAPPHSGRPGSQ were used as the immunogen for this GRIN3B antibody.

## Storage

Aliquot and store the GRIN3B antibody at -20oC.