


GRM7 Antibody (R34771)

Catalog No.	Formulation	Size
R34771-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2917
Applications	Western Blot : 0.5-2ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This GRM7 antibody is available for research use only.



Western blot testing of human cerebellum lysate with GRM7 antibody at 0.5ug/ml.
Predicted molecular weight: ~102kDa.

Description

Additional name(s) for this target protein: Glutamate receptor, metabotropic 7

Application Notes

Optimal dilution of the GRM7 antibody should be determined by the researcher.

Immunogen

Amino acids NCKLTISGSKKEDT were used as the immunogen for this GRM7 antibody.

Storage

Aliquot and store the GRM7 antibody at -20oC.