


Growth Hormone Antibody (pituitary) [clone GRHP1-1] (V7237)

Catalog No.	Formulation	Size
V7237-100UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide	100 ug
V7237-20UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide	20 ug
V7237SAF-100UG	1 mg/ml in 1X PBS; BSA free, sodium azide free	100 ug

[Bulk quote request](#)

Species Reactivity	Human
Format	Purified
Clonality	Monoclonal (mouse origin)
Isotype	Mouse IgG2b, kappa
Clone Name	GRHP1-1
Purity	Protein G affinity chromatography
Gene ID	2688
Localization	Cytoplasmic
Applications	Immunohistochemistry (FFPE) : 1-2ug/ml for 30 min at RT
Limitations	This Growth Hormone antibody is available for research use only.



IHC testing of FFPE human pituitary gland with Growth Hormone antibody (clone GRHP1-1). Required HIER: boil tissue sections in 10mM citrate buffer, pH 6, for 10-20 min.

Description

Pituitary growth hormone (GH) plays a crucial role in stimulating and controlling the growth, metabolism and differentiation of many mammalian cell types by modulating the synthesis of multiple mRNA species. These effects are mediated by the binding of GH to its membrane-bound receptor, GHR, and involve a phosphorylation cascade that results in the

modulation of numerous signaling pathways. GH is synthesized by acidophilic or somatotrophic cells of the anterior pituitary gland. Anti-GH is a useful marker in classification of pituitary tumors and the study of pituitary disease (acromegaly).

Application Notes

The concentration stated for each application is a general starting point. Variations in protocols, secondaries and substrates may require the Growth Hormone antibody to be titered up or down for optimal performance.

Immunogen

A human partial recombinant protein was used as the immunogen for this Growth Hormone antibody.

Storage

Store the Growth Hormone antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

References (1)