


HADH Antibody (R36362)


Catalog No.	Formulation	Size
R36362-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3033
Applications	Western Blot : 0.02-0.1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This HADH antibody is available for research use only.



Western blot testing of human kidney lysate with HADH antibody at 0.02ug/ml. Predicted molecular weight: ~33kDa.



Western blot of HEK293 lysate overexpressing HADH tested with HADH antibody (mock transfection in lane 1).

Description

Additional name(s) for this target protein: Hydroxyacyl-Coenzyme A dehydrogenase; L-3-hydroxyacyl-Coenzyme A dehydrogenase, short chain

Application Notes

Optimal dilution of the HADH antibody should be determined by the researcher.

Immunogen

Amino acids YERGDASKEDID were used as the immunogen for this HADH antibody.

Storage

Aliquot and store the HADH antibody at -20oC.