


Hap1 Antibody (R35211)

Catalog No.	Formulation	Size
R35211-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	15114
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This Hap1 antibody is available for research use only.



Western blot testing of mouse brain lysate with Hap1 antibody at 0.1ug/ml. Predicted molecular weight (mouse): ~67 kDa but can be observed at 67-100 kDa.

Description

Additional name(s) for this target protein: Huntingtin-associated protein 1

Application Notes

Optimal dilution of the Hap1 antibody should be determined by the researcher.

Immunogen

Amino acids QDAHSKRQKQK were used as the immunogen for this Hap1 antibody.

Storage

Aliquot and store the Hap1 antibody at -20oC.