


HCN3 Antibody (R35624)

Catalog No.	Formulation	Size
R35624-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Cow, Dog
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	57657
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This HCN3 antibody is available for research use only.



Western blot testing of rat kidney lysate with HCN3 antibody at 1ug/ml. Predicted molecular weight: ~86kDa.

Description

Additional name(s) for this target protein: Potassium/sodium hyperpolarization-activated cyclic nucleotide-gated channel 3

Application Notes

Optimal dilution of the HCN3 antibody should be determined by the researcher.

Immunogen

Amino acids RGVGRAPSTGAQ were used as the immunogen for this HCN3 antibody.

Storage

Aliquot and store the HCN3 antibody at -20oC.