


HHV8 ORF57 Antibody (F40036)


Catalog No.	Formulation	Size
F40036-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F40036-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	P88945
Applications	Western Blot : 1:1000
Limitations	This HHV8 ORF57 antibody is available for research use only.



HHV8 ORF57 antibody western blot analysis of over-expressed GFP-tagged ORF57 in HEK293T cell line.



HHV8 ORF57 antibody western blot analysis. ORF57 expressed from an inducible virally infected B-cell line, TREx BCBL 1-Rta.

Application Notes

Titration of the HHV8 ORF57 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 33-60 from the hlv8 protein was used as the immunogen for this HHV8 ORF57 antibody.

Storage

Aliquot the HHV8 ORF57 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.