


HIC2 Antibody (R35110)

Catalog No.	Formulation	Size
R35110-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Pig, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	23119
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This HIC2 antibody is available for research use only.



Western blot testing of HepG2 lysate with HIC2 antibody at 0.1ug/ml. Predicted molecular weight: ~66kDa.

Description

Additional name(s) for this target protein: Hypermethylated in cancer 2

Application Notes

Optimal dilution of the HIC2 antibody should be determined by the researcher.

Immunogen

Amino acids PQELPQAKGSDDE were used as the immunogen for this HIC2 antibody.

Storage

Aliquot and store the HIC2 antibody at -20oC.