


HIP2 Antibody / Huntingtin interacting protein 2 (R36108)


Catalog No.	Formulation	Size
R36108-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Dog, Pig, Cow
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3093
Applications	Western Blot : 1-3 ug/ml Immunohistochemistry (FFPE) : 1-2ug/ml ELISA (peptide) LOD : 1:2000
Limitations	This HIP2 antibody is available for research use only.



Western blot testing of Jurkat lysate with HIP2 antibody at 1.5ug/ml. Predicted molecular weight: ~22kDa.



IHC staining of FFPE human cortex with HIP2 antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: Huntingtin interacting protein 2, UBE2K

Application Notes

Optimal dilution of the HIP2 antibody should be determined by the researcher.

1. This HIP2 antibody is specific for isoforms 1, 2 and 3.

Immunogen

Amino acids SWDVETATELLSN were used as the immunogen for this HIP2 antibody.

Storage

Aliquot and store the HIP2 antibody at -20oC.