


## HIPPI Antibody (R33941)

| Catalog No.  | Formulation                                                                   | Size   |
|--------------|-------------------------------------------------------------------------------|--------|
| R33941-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

|                      |                                                         |
|----------------------|---------------------------------------------------------|
| Availability         | 1-3 business days                                       |
| Species Reactivity   | Human, Mouse                                            |
| Predicted Reactivity | Rat, Cow                                                |
| Format               | Antigen affinity purified                               |
| Clonality            | Polyclonal (goat origin)                                |
| Isotype              | Goat Ig                                                 |
| Purity               | Antigen affinity                                        |
| Gene ID              | 55081                                                   |
| Applications         | Western Blot : 1-3ug/ml<br>ELISA (peptide) LOD : 1:8000 |
| Limitations          | This HIPPI antibody is available for research use only. |



Western blot testing of mouse brain lysate with HIPPI antibody at 1ug/ml. Predicted molecular weight: ~49 kDa.

## Description

Additional name(s) for this target protein: HIP1 protein interactor; intraflagellar transport 57; IFT57

## Application Notes

Optimal dilution of the HIPPI antibody should be determined by the researcher.

## Immunogen

Amino acids EQPQEYDDPNATISN were used as the immunogen for this HIPPI antibody.

## Storage

Aliquot and store the HIPPI antibody at -20oC.