

HLA-DQ Antibody (MHC II) [clone SPV-L3] (V2581)


Catalog No.	Formulation	Size
V2581-100UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide	100 ug
V2581-20UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide	20 ug
V2581SAF-100UG	1 mg/ml in 1X PBS; BSA free, sodium azide free	100 ug




Citations (10)

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Clonality	Monoclonal (mouse origin)
Isotype	Mouse IgG2a, kappa
Clone Name	SPV-L3
Purity	Protein G affinity chromatography
UniProt	P01908, P01909, P01920
Localization	Cell surface
Applications	Flow Cytometry : 1-2ug/million cells Immunofluorescence : 1-2ug/ml
Limitations	This HLA-DQ antibody is available for research use only.



Flow cytometry testing of human Raji cells with HLA-DQ antibody (clone SPV-L3); Red=isotype control, Blue= HLA-DQ antibody.



SDS-PAGE analysis of purified, BSA-free HLA-DQ antibody (clone SPV-L3) as confirmation of integrity and purity.

Description

Recognizes a DQ antigen, which is a dimer of 60kDa. The class II molecule is a heterodimer consisting of an alpha (DQA) and a beta chain (DQB), both anchored in the membrane. It plays a central role in the immune system by presenting peptides derived from extracellular proteins. Class II molecules are expressed in antigen presenting cells (APC: B Lymphocytes, dendritic cells, macrophages). The alpha chain is approximately 33-35kDa. It is encoded by 5 exons; exon 1 encodes the leader peptide, exons 2 and 3 encode the two extracellular domains, and exon 4 encodes the transmembrane domain and the cytoplasmic tail. Within the DQ molecule both the alpha chain and the beta chain contain the polymorphisms specifying the peptide binding specificities, resulting in up to four different molecules. Typing for these polymorphisms is routinely done for bone marrow transplantation.

Application Notes

Optimal dilution of the HLA-DQ antibody should be determined by the researcher.

Immunogen

The T4-positive CTL clone HG-38 was used as the immunogen for the HLA-DQ antibody.

Storage

Store the HLA-DQ antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).