


HNRNPA1 Antibody / hnRNP A1 [clone ABEC-8] (RQ5286)

Catalog No.	Formulation	Size
RQ5286	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul

[Bulk quote request](#)

Availability	1-2 weeks
Species Reactivity	Human
Format	Purified
Clonality	Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	ABEC-8
Purity	Affinity purified
UniProt	P09651
Applications	Western Blot : 1:500-1:2000
Limitations	This HNRNPA1 antibody is available for research use only.



Western blot testing of human HepG2 cell lysate with HNRNPA1 antibody. Expected molecular weight: ~34 and/or 39 kDa.

Description

The HNRNPA1 gene encodes a member of a family of ubiquitously expressed heterogeneous nuclear ribonucleoproteins (hnRNPs), which are RNA-binding proteins that associate with pre-mRNAs in the nucleus and influence pre-mRNA processing, as well as other aspects of mRNA metabolism and transport. The protein encoded by this gene is one of the most abundant core proteins of hnRNP complexes and plays a key role in the regulation of alternative splicing. [RefSeq]

Application Notes

Optimal dilution of the HNRNPA1 antibody should be determined by the researcher.

Immunogen

A synthetic peptide specific to human hnRNP A1 / HNRNPA1 was used as the immunogen for the HNRNPA1 antibody.

Storage

Store the HNRNPA1 antibody at -20oC.