


HRASLS2 Antibody (R35247)

Catalog No.	Formulation	Size
R35247-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	54979
Applications	Immunohistochemistry (FFPE) : 5-10ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This HRASLS2 antibody is available for research use only.



IHC testing of FFPE human liver with HRASLS2 antibody at 5ug/ml. HIER: microwaved with pH9 Tris/EDTA buffer, HRP-staining.

Description

Additional name(s) for this target protein: HRAS-like suppressor 2

Application Notes

Optimal dilution of the HRASLS2 antibody should be determined by the researcher.

Immunogen

Amino acids QELPYSLTSDNC were used as the immunogen for this HRASLS2 antibody.

Storage

Aliquot and store the HRASLS2 antibody at -20oC.