


## HRH3 Antibody (R36104)

| Catalog No.  | Formulation                                                                   | Size   |
|--------------|-------------------------------------------------------------------------------|--------|
| R36104-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

|                      |                                                               |
|----------------------|---------------------------------------------------------------|
| Availability         | 1-3 business days                                             |
| Species Reactivity   | Human                                                         |
| Predicted Reactivity | Cow, Dog, Mouse, Rat                                          |
| Format               | Antigen affinity purified                                     |
| Clonality            | Polyclonal (goat origin)                                      |
| Isotype              | Goat Ig                                                       |
| Purity               | Antigen affinity                                              |
| UniProt              | Q9Y5N1                                                        |
| Gene ID              | 11255                                                         |
| Applications         | Western Blot : 0.1-0.3ug/ml<br>ELISA (peptide) LOD : 1:128000 |
| Limitations          | This HRH3 antibody is available for research use only.        |



Western blot testing of A549 lysate with HRH3 antibody at 0.1ug/ml. Predicted molecular weight: 36~49kDa (isoforms 1-7).

## Description

Additional name(s) for this target protein: Histamine H3 receptor

## Application Notes

Optimal dilution of the HRH3 antibody should be determined by the researcher.

## Immunogen

Amino acids SVTRAVSYRAQQGDT were used as the immunogen for this HRH3 antibody. This sequence is common to isoforms 1-7 (UniProt Q9Y5N1).

## Storage

Aliquot and store the HRH3 antibody at -20oC.