


HRS Antibody (R34908)

Catalog No.	Formulation	Size
R34908-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	9146
Applications	Immunohistochemistry (FFPE) : 1-2ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This HRS antibody is available for research use only.



IHC testing of FFPE human testis with HRS antibody at 2.5ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: Hepatocyte growth factor-regulated tyrosine kinase substrate; HGS

Application Notes

Optimal dilution of the HRS antibody should be determined by the researcher.

Immunogen

Amino acids NRNYWEKKQEEARK were used as the immunogen for this HRS antibody.

Storage

Aliquot and store the HRS antibody at -20oC.