


HSF1 Antibody (F48065)

Catalog No.	Formulation	Size
F48065-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F48065-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Bovine
Format	Purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Purified
UniProt	Q00613
Applications	Immunofluorescence : 1:10-1:50
Limitations	This HSF1 antibody is available for research use only.



Confocal immunofluorescent analysis of HSF1 antibody with HeLa cells followed by Alexa Fluor 488-conjugated goat anti-rabbit IgG (green). DAPI was used as a nuclear counterstain (blue).

Description

Heat shock transcription factor 1 (HSF1) mediates the induction of heat shock protein gene expression in cells exposed to elevated temperature and other stress conditions. In response to stress, HSF1 acquires DNA-binding ability and localizes to nuclear stress granules. SUMO modification of HSF1 converts HSF1 to the DNA-binding form. HSF1 colocalizes with SUMO-1 in nuclear stress granules, which is prevented by mutation of the HSF1 lysine targeted for sumoylation.

Application Notes

Titration of the HSF1 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 278-309 from the human protein (sumoylation site) was used as the immunogen for this HSF1 antibody.

Storage

Aliquot the HSF1 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.