


## HSP27 Antibody [clone HII-8] (RQ5080)

Catalog No.	Formulation	Size
RQ5080	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul

[Bulk quote request](#)

Availability	1-2 weeks
Species Reactivity	Human
Format	Purified
Clonality	Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	HII-8
Purity	Affinity purified
UniProt	P04792
Applications	Western Blot : 1:500-1:1000
Limitations	This HSP27 antibody is available for research use only.



Western blot testing of human HeLa cell lysate with HSP27 antibody. Expected molecular weight: 23-27 kDa.

### Description

Heat shock protein beta-1, also called HSP27, is a small heat shock protein which functions as a molecular chaperone probably maintaining denatured proteins in a folding-competent state. Plays a role in stress resistance and actin organization. Through its molecular chaperone activity may regulate numerous biological processes including the phosphorylation and the axonal transport of neurofilament proteins. [UniProt]

### Application Notes

Optimal dilution of the HSP27 antibody should be determined by the researcher.

## **Immunogen**

A synthetic peptide specific to human HSP27 / HSPB1 was used as the immunogen for the HSP27 antibody.

## **Storage**

Store the HSP27 antibody at -20oC.