


HSP90AA1 Antibody / HSP90 alpha [clone BAD-8] (RQ5192)

Catalog No.	Formulation	Size
RQ5192	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul

[Bulk quote request](#)

Availability	1-2 weeks
Species Reactivity	Human
Format	Purified
Clonality	Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	BAD-8
Purity	Affinity purified
UniProt	P07900
Applications	Western Blot : 1:500-1:2000
Limitations	This HSP90AA1 antibody is available for research use only.



Western blot testing of human Jurkat cell lysate with HSP90AA1 antibody. Predicted molecular weight: 86-90 kDa.

Description

The protein encoded by the HSP90AA1 gene is an inducible molecular chaperone that functions as a homodimer. The encoded protein aids in the proper folding of specific target proteins by use of an ATPase activity that is modulated by co-chaperones. [RefSeq]

Application Notes

Optimal dilution of the HSP90AA1 antibody should be determined by the researcher.

Immunogen

A synthetic peptide specific to human HSP90 alpha / HSP90AA1 was used as the immunogen for the HSP90AA1 antibody.

Storage

Store the HSP90AA1 antibody at -20oC.