


HSPA4L Antibody / Heat shock 70 kDa protein 4L (R36492)

Catalog No.	Formulation	Size
R36492-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
UniProt	O95757
Localization	Cytoplasmic, nuclear
Applications	IHC (FFPE) : 5-10ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This HSPA4L antibody is available for research use only.



IHC testing of FFPE human testis tissue with HSPA4L antibody at 5ug/ml. HIER: steamed antigen retrieval with pH6 citrate buffer; AP-staining.

Description

Alternate name(s) for this target: Heat shock 70kDa protein 4-like; APG-1; Osp94; heat shock 70 kDa protein 4-like protein; heat shock protein (hsp110 family)

Application Notes

Optimal dilution of the HSPA4L antibody should be determined by the researcher.

Immunogen

Amino acids DSSQHTKSSGEMEVD were used as the immunogen for this HSPA4L antibody.

Storage

Aliquot and store the HSPA4L antibody at -20oC.