


HSPE1 Antibody (R35007)

Catalog No.	Formulation	Size
R35007-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Pig
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3336
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This HSPE1 antibody is available for research use only.



Western blot testing of A) HeLa and B) HepG2 lysates with HSPE1 antibody at 0.3ug/ml. Expected molecular weight: ~10kDa.

Description

Additional name(s) for this target protein: Heat shock protein 10, HSP10

Application Notes

Optimal dilution of the HSPE1 antibody should be determined by the researcher.

Immunogen

Amino acids PEYGGTKVVLDDK were used as the immunogen for this HSPE1 antibody.

Storage

Aliquot and store the HSPE1 antibody at -20oC.