


HTR2C Antibody (R35501)

Catalog No.	Formulation	Size
R35501-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat, Dog, Cow
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3358
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This HTR2C antibody is available for research use only.



Western blot testing of EBV immortalized lymphoblastoid lysate with HTR2C antibody at 1ug/ml. Predicted molecular weight: ~52 kDa, but may be observed at higher molecular weights due to glycosylation.

Description

Additional name(s) for this target protein: 5-hydroxytryptamine (serotonin) receptor 2C,

Application Notes

Optimal dilution of the HTR2C antibody should be determined by the researcher.

Immunogen

Amino acids QVENLELPVN were used as the immunogen for this HTR2C antibody.

Storage

Aliquot and store the HTR2C antibody at -20oC.