


ICAM2 Antibody (RQ4133)


Catalog No.	Formulation	Size
RQ4133	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human, Mouse
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Antigen affinity purified
Buffer	Lyophilized from 1X PBS with 2% Trehalose
UniProt	P13598
Localization	Cell membrane
Applications	Western Blot : 0.5-1ug/ml Flow Cytometry : 1-3ug/million cells Direct ELISA : 0.1-0.5ug/ml
Limitations	This ICAM2 antibody is available for research use only.



Flow cytometry testing of fixed human HEL cells with ICAM2 antibody at 1ug/million cells (blocked with goat sera); Red=cells alone, Green=isotype control, Blue= ICAM2 antibody.



Western blot testing of 1) human Jurkat, 2) human HEL, 3) human K562 and 4) mouse spleen tissue lysate with ICAM2 antibody at 0.5ug/ml. Expected molecular weight: 31-60 kDa depending on glycosylation level.

