


IDE Antibody Biotin Conjugate (R35417BTN)

Catalog No.	Formulation	Size
R35417BTN	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog
Format	Biotin Conjugate
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3416
Applications	Western Blot : 0.3-1ug/ml Immunohistochemistry (FFPE) : suitable ELISA (peptide) LOD : 1:32000
Limitations	This IDE antibody is available for research use only.



Western blot testing of K562 cell lysate with biotinylated IDE antibody at 0.3ug/ml.
Predicted molecular weight: ~118kDa.

Description

Additional name(s) for this target protein: Insulin-degrading enzyme, Insulysin, Insulinase

Application Notes

Optimal dilution of the IDE antibody should be determined by the researcher.

Immunogen

Amino acids QQYNFDRDNTE were used as the immunogen for this IDE antibody.

Storage

Aliquot and store the IDE antibody at -20oC.