

IER3 Antibody / Immediate early response gene 3 (FY12932)

Catalog No.	Formulation	Size
FY12932	Adding 0.2 ml of distilled water will yield a concentration of 500 ug/ml	100 ug

[Bulk quote request](#)

Availability	1-2 days
Species Reactivity	Human
Format	Lyophilized
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Immunogen affinity purified
Buffer	Each vial contains 4 mg Trehalose, 0.9 mg NaCl, 0.2 mg Na ₂ HPO ₄ .
UniProt	P46695
Applications	Western Blot : 0.25-0.5ug/ml Immunohistochemistry : 2-5ug/ml Flow Cytometry : 1-3ug/million cells
Limitations	This IER3 antibody is available for research use only.

Description

IER3 antibody detects Immediate early response gene 3 protein, a stress-inducible regulator of apoptosis and cell survival pathways. The UniProt recommended name is Immediate early response 3 (IER3), also known as IEX-1, DIF-2, gly96, and radiation-inducible protein. IER3 acts as a rapid responder to stress, mitogenic, and inflammatory signals, mediating downstream responses to NF-kappaB, p53, and MAPK pathways. Its expression is transiently induced by stimuli such as UV radiation, cytokines (TNF-alpha, IL-1), and growth factors, reflecting its role in adaptive cell signaling.

Functionally, IER3 antibody targets a multifunctional protein that modulates apoptosis by interacting with BCL2 family members and regulating mitochondrial membrane potential. It can act as either pro- or anti-apoptotic depending on cell context and stimulus. IER3 has been implicated in cell cycle control, stress tolerance, and immune regulation, with overexpression protecting cells from TNF-induced apoptosis via suppression of caspase-3 activation. Conversely, downregulation promotes sensitivity to radiation and chemotherapeutic stress, linking IER3 to cancer therapy resistance mechanisms.

IER3 localizes mainly to the cytoplasm but can also translocate to the nucleus in certain stress conditions. The IER3 antibody is valuable for monitoring expression patterns across tissues, including spleen, thymus, and epithelial cells, and

in cancers such as leukemia, breast, and colon carcinoma. The IER3 gene is located on chromosome 6p21.33, close to the major histocompatibility complex (MHC) region, encoding a 156-amino acid protein. Its promoter contains NF-kappaB binding sites, allowing rapid induction following stress or infection. IER3 acts as a negative feedback regulator in immune responses, balancing cytokine signaling and cellular resilience.

IER3 also regulates ERK and AKT phosphorylation states, controlling cellular energy metabolism and proliferation. Its modulation of mitochondrial ATP synthase interaction contributes to ROS balance. NSJ Bioreagents provides IER3 antibody reagents optimized for detecting both native and induced protein forms across species for use in western blotting, immunocytochemistry, and flow cytometry applications.

Through its unique duality in regulating cell survival, IER3 remains a critical molecular switch in stress adaptation, immune tolerance, and oncogenesis.

Application Notes

Optimal dilution of the IER3 antibody should be determined by the researcher.

Immunogen

A synthetic peptide corresponding to a sequence at the C-terminus of human IER3 was used as the immunogen for the IER3 antibody.

Storage

After reconstitution, the IER3 antibody can be stored for up to one month at 4oC. For long-term, aliquot and store at -20oC. Avoid repeated freezing and thawing.