


IFI204 Antibody (N1138)

Catalog No.	Formulation	Size
N1138-100UG	0.5 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide	100 ug
N1138-25UG	0.5 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide	25 ug

[Bulk quote request](#)

Species Reactivity	Mouse
Format	Purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Protein A affinity chromatography purified antibody
Buffer	1X PBS, pH 7.4
Gene ID	15951
Applications	Western Blot : 3-6ug/ml
Limitations	This IFI204 antibody is available for research use only.



Western blot testing of mouse 1) kidney and 2) EL4 lysate with IFI204 antibody at 3ug/ml.

Description

Part of the innate immune response to infection is the production of type I IFNs (interferons). The murine Ifi200 family of proteins are activated by IFNs. IFI204 (Interferon-activable protein 204), or p204, is a regulator of bone production and the differentiation of macrophages and skeletal and cardiac muscle cells.

Application Notes

Provided assay concentrations are suggestions only, antibody titration may be required for optimal results.

Immunogen

A recombinant protein fragment from the N-terminal region of mouse IFI204 was used as the immunogen for this

antibody.

Storage

Aliquot and store frozen at -20oC or colder to avoid repeated freeze-thaw cycles.

Alternate Names

interferon activated gene 204, IFI16 (human)