


Ihh Antibody (F40064)

Catalog No.	Formulation	Size
F40064-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F40064-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q80XI9
Applications	IHC (Paraffin) : 1:10-1:50 Western Blot : 1:1000
Limitations	This Ihh antibody is available for research use only.



IHC analysis of FFPE human lung carcinoma tissue stained with Ihh antibody



Ihh antibody western blot analysis in mouse liver tissue lysate

Description

Ihh (Indian hedgehog) is a critical and possibly direct regulator of joint development. In its absence, distribution and function of Gdf5-expressing interzone-associated cells are abnormal.

Application Notes

Titration of the Ihh antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 1-30 from the mouse protein was used as the immunogen for this Ihh antibody.

Storage

Aliquot the Ihh antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.