


Ihh Antibody (F40065)

Catalog No.	Formulation	Size
F40065-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F40065-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Chicken, Zebrafish, Xenopus
Format	Purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Purified
UniProt	Q80XI9
Applications	Western Blot : 1:1000
Limitations	This Ihh antibody is available for research use only.



Description

Ihh (Indian hedgehog) is a critical and possibly direct regulator of joint development. In its absence, distribution and function of Gdf5-expressing interzone-associated cells are abnormal.

Application Notes

Titration of the Ihh antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 169-200 from the mouse protein was used as the immunogen for this Ihh antibody.

Storage

Aliquot the Ihh antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.