



IkB alpha Antibody (F50935)

Catalog No.	Formulation	Size
F50935-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F50935-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml


[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Pig
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	P25963
Applications	Western Blot : 1:1000 IHC (Paraffin) : 1:50-1:100 Flow Cytometry : 1:10-1:50 Immunofluorescence : 1:10-1:50
Limitations	This IkB alpha antibody is available for research use only.






IHC analysis of FFPE human hepatocarcinoma with IκB alpha antibody



Flow cytometric analysis of HeLa cells using IκB alpha antibody (bottom histogram) compared to a negative control (top histogram). FITC-conjugated goat-anti-rabbit secondary Ab was used for the analysis.



Confocal immunofluorescent analysis of IκB alpha antibody with ZR-75-1 cells followed by Alexa Fluor 488-conjugated goat anti-rabbit IgG (green). Actin filaments have been labeled with Alexa Fluor 555 Phalloidin (red).

Description

NFKB1 or NFKB2 is bound to REL, RELA, or RELB to form the NFKB complex. The NFKB complex is inhibited by I-κappa-B proteins (NFKBIA or NFKBIB), which inactivate NF-κappa-B by trapping it in the cytoplasm. Phosphorylation of serine residues on the I-κappa-B proteins by kinases (IKBKA, IKBKB) marks them for destruction via the ubiquitination pathway, thereby allowing activation of the NF-κappa-B complex. Activated NFKB complex translocates into the nucleus and binds DNA at kappa-B-binding motifs such as 5-prime GGGRNYYCC 3-prime or 5-prime HGGARNYYCC 3-prime (where H is A, C, or T; R is an A or G purine; and Y is a C or T pyrimidine).

Application Notes

Titration of the IκB alpha antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 12-41 from the human protein was used as the immunogen for this IκB alpha antibody.

Storage

Aliquot the IκB alpha antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.

