


IL-1B Antibody (F51185)


Catalog No.	Formulation	Size
F51185-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F51185-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)


Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Primate, Rabbit
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	P01584
Applications	Western Blot : 1:1000 IHC (Paraffin) : 1:10-1:50 Flow Cytometry : 1:10-1:50
Limitations	This IL-1B antibody is available for research use only.



Flow cytometric analysis of NCI-H292 cells using IL-1B antibody (bottom histogram) compared to a negative control (top histogram). FITC-conjugated goat-anti-rabbit secondary Ab was used for the analysis.



IHC analysis of FFPE human lung carcinoma with IL-1B antibody



95
72
55
43
34
26

Western blot analysis of IL-1B antibody and NCI-H460 lysate. Predicted molecular weight: 31 kDa (precursor) and 17 kDa (mature).

Description

IL1B is produced by activated macrophages, IL-1 stimulates thymocyte proliferation by inducing IL-2 release, B-cell maturation and proliferation, and fibroblast growth factor activity. IL-1 proteins are involved in the inflammatory response, being identified as endogenous pyrogens, and are reported to stimulate the release of prostaglandin and collagenase from synovial cells. [UniProt]

Application Notes

Titration of the IL-1B antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 148-174 from the human protein was used as the immunogen for this IL-1B antibody.

Storage

Aliquot the IL-1B antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.