


IL-2 Antibody (R30538)

Catalog No.	Formulation	Size
R30538	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Mouse
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Antigen affinity
Buffer	Lyophilized from 1X PBS with 2.5% BSA and 0.025% sodium azide/thimerosal
UniProt	P04351
Applications	Western Blot : 0.5-1ug/ml
Limitations	This IL-2 antibody is available for research use only.



Western blot testing of IL-2 antibody and Lane 1: recombinant mouse protein 10ng; 2: 5ng; 3: 2.5ng

Description

Interleukin-2, formerly referred to as T-cell growth factor, is a powerfully immunoregulatory lymphokine that is produced by lectin-or antigen-activated T cells. The protein can act as a growth hormone for both B and T lymphocytes. It is useful in the study of the molecular nature of T-cell differentiation and, like interferons, augments natural killer cell activity.

Application Notes

The stated application concentrations are suggested starting amounts. Titration of the IL-2 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

An amino acid sequence from the C-terminus of mouse IL-2 (DFLRRWIAFCQSIISTSPQ) was used as the immunogen for this IL-2 antibody.

Storage

After reconstitution, the IL-2 antibody can be stored for up to one month at 4°C. For long-term, aliquot and store at -20°C. Avoid repeated freezing and thawing.