


IL-22 Antibody (F43455)

Catalog No.	Formulation	Size
F43455-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F43455-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q9GZX6
Localization	Secreted
Applications	Western Blot : 1:1000
Limitations	This IL-22 antibody is available for research use only.



Description

Cytokine that contributes to the inflammatory response in vivo. [UniProt]

Application Notes

Titration of the IL-22 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 44-73 from the human protein was used as the immunogen for this IL-22 antibody.

Storage

Aliquot the IL-22 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.