


IL-22 Antibody (R34201)

Catalog No.	Formulation	Size
R34201-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	53342
Localization	Secreted
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This IL-22 antibody is available for research use only.



Western blot testing of human frontal cortex lysate with IL-22 antibody at 1ug/ml.
Observed molecular weight: 19~25 kDa depending on glycosylation level.

Description

Additional name(s) for this target protein: Interleukin-17D, interleukin-27

Application Notes

Optimal dilution of the IL-22 antibody should be determined by the researcher.

Immunogen

Amino acids HHTLQLGPREQ were used as the immunogen for this IL-22 antibody.

Storage

Aliquot and store the IL-22 antibody at -20oC.