


Il12b Antibody (R33713)

Catalog No.	Formulation	Size
R33713-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	16160
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This Il12b antibody is available for research use only.



Western blot testing of 1) mouse skin and 2) rat skin lysate with Il12b antibody at 0.5ug/ml. Predicted molecular weight: 37~40 kDa.

Description

Additional name(s) for this target protein: Interleukin 12B

Application Notes

Optimal dilution of the Il12b antibody should be determined by the researcher.

Immunogen

Amino acids EARQQNKYENYS were used as the immunogen for this Il12b antibody.

Storage

Aliquot and store the IL12b antibody at -20oC.