


IL12B Antibody (R35275)

Catalog No.	Formulation	Size
R35275-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3593
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This IL12B antibody is available for research use only.



Western blot testing of human 1) liver, 2) skin and 3) tonsil lysate with IL12B antibody at 0.3ug/ml. Predicted molecular weight: 37~40 kDa.

Description

Additional name(s) for this target protein: Interleukin 12B (natural killer cell stimulatory factor 2, cytotoxic lymphocyte maturation factor 2, p40); IL12/p40

Application Notes

Optimal dilution of the IL12B antibody should be determined by the researcher.

Immunogen

Amino acids QGKSKREKKDRVFTD were used as the immunogen for this IL12B antibody.

Storage

Aliquot and store the IL12B antibody at -20oC.