


IL1RN Antibody / IL1RA / Interleukin-1 Receptor Antagonist (F54851)


Catalog No.	Formulation	Size
F54851-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F54851-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

Bulk quote request


Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Purified
UniProt	P18510
Localization	Cytoplasmic
Applications	Flow Cytometry : 1:10-1:50 (1x10 ⁶ cells) Immunofluorescence : 1:10-1:50 Immunohistochemistry (FFPE) : 1:10-1:50 Western Blot : 1:500-1:1000
Limitations	This IL1RN antibody is available for research use only.




IHC testing of FFPE human colon carcinoma tissue with IL1RN antibody. HIER: steam section in pH6 citrate buffer for 20 min and allow to cool prior to staining.



Immunofluorescent staining of human WiDr cells with IL1RN antibody (green) and DAPI nuclear stain (blue).



Western blot testing of 1) non-transfected and 2) transfected 293 cell lysate with IL1RN antibody.



Flow cytometry testing of human HEK293 cells with IL1RN antibody; Blue=isotype control, Green= IL1RN antibody.

Description

The protein is a member of the interleukin 1 cytokine family. This protein inhibits the activities of interleukin 1, alpha (IL1A) and interleukin 1, beta (IL1B), and modulates a variety of interleukin 1 related immune and inflammatory responses. A polymorphism of its gene is reported to be associated with increased risk of osteoporotic fractures and gastric cancer.

Application Notes

The stated application concentrations are suggested starting points. Titration of the IL1RN antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 100-129 from the human protein was used as the immunogen for the IL1RN antibody.

Storage

Aliquot the IL1RN antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.